



Reso-Strip[®]


Sterile, self-adhesive
skin-closure tapes made of
polyamide for better cosmetic results

 **RESORBA[®]**
REPARIEREN UND REGENERIEREN

LATEX-FREE

Applications:

- ▶ Supports and effectively relieves the tension on subcuticular, interrupted or staple sutures
- ▶ To relieve the tension on fresh wounds where the sutures have been removed early to achieve better cosmetic results
- ▶ As a form of wound closure beneath support dressings
- ▶ To approximate the wound margins of wounds healing by secondary intention
- ▶ No visible sutures when placed on wounds closed with subcutaneous sutures
- ▶ For the atraumatic closure of wounds in regions where a good cosmetic result is of particular importance



Figures showing original sizes

Product qualities:

- ▶ Rounded corners guarantee excellent adherence
- ▶ Fast and easy handling
- ▶ Hypoallergenic acrylate adhesive, dermatologically tested
- ▶ Painless application and removal without the need for local anaesthetic
- ▶ Encourages the healing process through its high level of air permeability
- ▶ The dispensing package can be inserted without problem into the usual suture rack
- ▶ Its skin colour renders it inconspicuous on the skin

Product range:

| | REF | Dimensions | Contents/sachet | Sachets/pack |
|--------------------|------|-------------|-----------------|--------------|
| Reso-Strip® | 200 | 3 x 76 mm | 6 | 36 |
| | 201 | 6 x 76 mm | 6 | 36 |
| | 202 | 6 x 38 mm | 6 | 36 |
| | 203 | 6 x 101 mm | 10 | 12 |
| | 204 | 12 x 101 mm | 6 | 12 |
| | 205 | 25 x 125 mm | 4 | 12 |
| | 206 | 12 x 50 mm | 4 | 12 |
| Praxi-Pack | 1201 | 6 x 76 mm | 3 | 36 |
| | 1202 | 6 x 38 mm | 3 | 36 |
| | 1204 | 12 x 101 mm | 3 | 12 |

Resorba Wundversorgung GmbH + Co. KG
Am Flachmoor 16
D-90475 Nürnberg

Fon: +49 (0) 91 28 - 9 11 5 - 0
Fax: +49 (0) 91 28 - 9 11 5 - 91
E-Mail: infomail@resorba.com
www.resorba.com

CE 1275

08/2009 P1077

RESORBA®
REPARIEREN UND REGENERIEREN