











# CAPROLON RESORBA® P (LA/CL)

Sterile Nadel-Faden-Kombination









Needle - suture - combination, sterile

| Nadel / needle   | Stärke/gauge |     | Länge/length<br>m | Farbe/colour   | Best.-Nr.<br>code no. | Pckg.-Inhalt<br>box content<br>in Dtzd./dozen |
|--|--------------|-----|-------------------|----------------|-----------------------|---|
|  | metric       | USP |                   |                |                       |   |
| <b>5/8-Kreis Rundkörper = FR</b><br>5/8 circle, round bodied = FR                                  |              |     |                   |                |                       |   |
| ● 2 x FR 27<br>   | 3            | 2-0 | 0,70              | violett/violet | PM 2011               | 2   |
| <b>1/2-Kreis Rundkörper = HR, HRS</b><br>1/2 circle, round bodied = HR, HRS                        |              |     |                   |                |                       |   |
| ● HR 12<br>       | 0.7          | 6-0 | 0,70              | violett/violet | PM 2021               | 2   |
|  | 1            | 5-0 | 0,70              | violett/violet | PM 2022               | 2   |
|  | 1.5          | 4-0 | 0,70              | violett/violet | PM 2023               | 2   |
| ● 2 x HR 12<br>  | 1            | 5-0 | 0,70              | violett/violet | PM 2024               | 2   |
|  | 1.5          | 4-0 | 0,70              | violett/violet | PM 20222              | 2   |
| ● HRS 12<br>    | 2            | 3-0 | 0,70              | violett/violet | PM 2031               | 2   |
|  | 3            | 2-0 | 0,70              | violett/violet | PM 2032               | 2   |
|  | 3.5          | 0   | 0,70              | violett/violet | PM 2033               | 2   |
| ● HR 17<br>     | 0.7          | 6-0 | 0,70              | violett/violet | PM 2025               | 2   |
|  | 1            | 5-0 | 0,70              | violett/violet | PM 2026               | 2   |
|  | 1.5          | 4-0 | 0,70              | violett/violet | PM 2027               | 2   |
|  | 2            | 3-0 | 0,70              | violett/violet | PM 2028               | 2   |
| ● 2 x HR 17<br> | 1            | 5-0 | 0,70              | violett/violet | PM 2029               | 2   |
|  | 1.5          | 4-0 | 0,70              | violett/violet | PM 20210              | 2   |
| ● HR 22<br>     | 1.5          | 4-0 | 0,70              | violett/violet | PM 20211              | 2   |
|  | 2            | 3-0 | 0,70              | violett/violet | PM 20212              | 2   |
| ● HRS 22<br>    | 3            | 2-0 | 0,70              | violett/violet | PM 2034               | 2   |
|  | 3.5          | 0   | 0,70              | violett/violet | PM 2035               | 2   |
| ● HR 27<br>     | 1.5          | 4-0 | 0,70              | violett/violet | PM 20213              | 2   |
|  | 2            | 3-0 | 0,70              | violett/violet | PM 20214              | 2   |
|  | 3            | 2-0 | 0,70              | violett/violet | PM 20215              | 2   |
|  | 2            | 3-0 | 6 x 0,70          | violett/violet | PM 20216              | 1   |
| ● 2 x HR 27<br> | 1.5          | 4-0 | 0,70              | violett/violet | PM 20217              | 2   |
|  | 2            | 3-0 | 0,70              | violett/violet | PM 20218              | 2   |
|  | 3            | 2-0 | 0,70              | violett/violet | PM 20219              | 2   |

# CAPROLON RESORBA® P (LA/CL)

Sterile Nadel-Faden-Kombination










Needle - suture - combination, sterile

| Nadel / needle  | Stärke/gauge |     | Länge/length<br>m | Farbe/colour   | Best.-Nr.<br>code no. | Pckg.-Inhalt<br>box content<br>in Dtzd./dozen |
|---|--------------|-----|-------------------|----------------|-----------------------|---|
|   | metric       | USP |                   |                |                       |   |
| ● HRS 27<br>                     | 3            | 2-0 | 0,70              | violett/violet | PM 2036               | 2   |
|   | 3.5          | 0   | 0,70              | violett/violet | PM 2037               | 2   |
|   | 4            | 1   | 0,70              | violett/violet | PM 2038               | 2   |
| ● HR 35<br>                      | 2            | 3-0 | 0,70              | violett/violet | PM 20220              | 2   |
|   | 3            | 2-0 | 0,70              | violett/violet | PM 20221              | 2   |
| ● HRS 35<br>                     | 3            | 2-0 | 0,70              | violett/violet | PM 2039               | 2   |
|   | 3.5          | 0   | 0,70              | violett/violet | PM 20310              | 2   |
|   | 4            | 1   | 0,70              | violett/violet | PM 20311              | 2   |
| ● HRS 40<br>                    | 3            | 2-0 | 0,90              | violett/violet | PM 20312              | 2   |
|   | 3.5          | 0   | 0,90              | violett/violet | PM 20313              | 2   |
|   | 4            | 1   | 0,90              | violett/violet | PM 20314              | 2   |
| ● HRS 40 / Schlinge / loop<br> | 3.5          | 0   | 1,50              | violett/violet | PM 20318              | 2   |
|   | 4            | 1   | 1,50              | violett/violet | PM 20319              | 2   |
| ● HRS 48<br>                   | 3.5          | 0   | 0,90              | violett/violet | PM 20315              | 2   |
|   | 4            | 1   | 0,90              | violett/violet | PM 20316              | 2   |
|   | 5            | 2   | 0,90              | violett/violet | PM 20317              | 2   |
| ● 2 x HRS 48<br>               | 4            | 1   | 0,90              | violett/violet | PM 20321              | 2   |
| ● HRS 48 / Schlinge / loop<br> | 4            | 1   | 1,80              | violett/violet | PM 20320              | 2   |

# CAPROLON RESORBA® P (LA/CL)

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

| Nadel / needle   | Stärke/gauge |     | Länge/length<br>m | Farbe/colour   | Best.-Nr.<br>code no. | Pckg.-Inhalt<br>box content<br>in Dtzd./dozen |
|--|--------------|-----|-------------------|----------------|-----------------------|---|
|  | metric       | USP |                   |                |                       |   |
| <b>1/2-Kreis schneidende Rundkörper = HRT</b><br><i>1/2-circle round bodied, trocar point = HRT</i>                      |              |     |                   |                |                       |   |
| <b>▼ HRT 18</b><br>                     | 1            | 5-0 | 0,70              | violett/violet | PM 2061               | 2   |
|  | 1.5          | 4-0 | 0,70              | violett/violet | PM 2062               | 2   |
|  | 2            | 3-0 | 0,70              | violett/violet | PM 2063               | 2   |
|  | 3            | 2-0 | 0,70              | violett/violet | PM 2064               | 2   |
| <b>▼ HRT 22</b><br>                     | 1.5          | 4-0 | 0,70              | violett/violet | PM 2065               | 2   |
|  | 2            | 3-0 | 0,70              | violett/violet | PM 2066               | 2   |
| <b>▼ HRT 25</b><br>                     | 1.5          | 4-0 | 0,70              | violett/violet | PM 2069               | 2   |
|  | 2            | 3-0 | 0,70              | violett/violet | PM 2067               | 2   |
|  | 3            | 2-0 | 0,70              | violett/violet | PM 2068               | 2   |
| <b>▼ HRT 50 / Schlinge / loop</b><br> | 4            | 1   | 1,80              | violett/violet | PM 20610              | 2   |
| <b>3/8-Kreis Rundkörper = DR</b><br><i>3/8 circle, round bodied = DR</i>   |              |     |                   |                |                       |   |
| <b>● 2 x DR 10</b><br>                | 0.5          | 7-0 | 0,70              | violett/violet | PM 2071               | 2   |
|  | 0.7          | 6-0 | 0,70              | violett/violet | PM 2072               | 2   |
| <b>● DR 12</b><br>                    | 0.7          | 6-0 | 0,70              | violett/violet | PM 2073               | 2   |
|  | 1            | 5-0 | 0,70              | violett/violet | PM 2074               | 2   |
| <b>● 2 x DR 12</b><br>                | 0.7          | 6-0 | 0,70              | violett/violet | PM 2075               | 2   |
|  | 1            | 5-0 | 0,70              | violett/violet | PM 2076               | 2   |
| <b>Gerade Rundkörper = GR</b><br><i>straight, round bodied = GR</i>  |              |     |                   |                |                       |   |
| <b>● 2 x GR 20</b><br>                | 1.5          | 4-0 | 0,35              | violett/violet | PM 2091               | 2   |
| <b>● 2 x GR 65</b><br>                | 3.5          | 0   | 0,70              | violett/violet | PM 2092               | 2   |
|  | 4            | 1   | 0,70              | violett/violet | PM 2093               | 2   |

# CAPROLON RESORBA® P (LA/CL)

Sterile Nadel-Faden-Kombination









Needle - suture - combination, sterile

| Nadel / needle  | Stärke/gauge |     | Länge/length<br>m | Farbe/colour         | Best.-Nr.<br>code no. | Pckg.-Inhalt<br>box content<br>in Dtzd./dozen |
|---|--------------|-----|-------------------|----------------------|-----------------------|---|
|   | metric       | USP |                   |                      |                       |   |
| <b>1/2-Kreis schneidende = HS</b><br><i>1/2 circle, reverse cutting = HS</i>  |              |     |                   |                      |                       |   |
| ▼ <b>HS 27</b><br>                                     | 3            | 2-0 | 0,70              | violett/violet       | PM 2111               | 2   |
|   | 3.5          | 0   | 0,70              | violett/violet       | PM 2112               | 2   |
|   | 4            | 1   | 0,70              | violett/violet       | PM 2113               | 2   |
| ▼ <b>HS 35</b><br>                                     | 3.5          | 0   | 0,70              | violett/violet       | PM 2114               | 2   |
|   | 4            | 1   | 0,70              | violett/violet       | PM 2115               | 2   |
| ▼ <b>HS 40</b><br>                                     | 3.5          | 0   | 0,90              | violett/violet       | PM 2116               | 2   |
| <b>3/8-Kreis schneidende = DS,<br/>mit PREMIUM-Schliff = DSM</b><br><i>3/8 circle, reverse cutting = DS,<br/>PREMIUM cut = DSM</i>      |              |     |                   |                      |                       |   |
| ▼ <b>DSM 11</b><br><br>PREMIUM-Schliff / PREMIUM cut | 0.7          | 6-0 | 0,45              | violett/violet       | PM 2151               | 2   |
|   | 0.7          | 6-0 | 0,45              | ungefärbt/<br>undyed | PM 2152               | 2   |
|   | 1            | 5-0 | 0,45              | ungefärbt/<br>undyed | PM 2153               | 2   |
| ▽ <b>DSM 13</b><br><br>PREMIUM-Schliff / PREMIUM cut | 1            | 5-0 | 0,45              | violett/violet       | PM 2154               | 2   |
|   | 1.5          | 4-0 | 0,45              | violett/violet       | PM 2155               | 2   |
|   | 1            | 5-0 | 0,45              | ungefärbt/<br>undyed | PM 2156               | 2   |
|   | 1.5          | 4-0 | 0,45              | ungefärbt/<br>undyed | PM 2157               | 2   |
| ▽ <b>DS 18</b><br>                                   | 1            | 5-0 | 0,45              | ungefärbt/<br>undyed | PM 2146               | 2   |
|   | 1.5          | 4-0 | 0,45              | ungefärbt/<br>undyed | PM 2141               | 2   |
|   | 2            | 3-0 | 0,45              | ungefärbt/<br>undyed | PM 2142               | 2   |
| ▽ <b>DSM 18</b><br><br>PREMIUM-Schliff / PREMIUM cut | 1.5          | 4-0 | 0,45              | violett/violet       | PM 21510              | 2   |
|   | 2            | 3-0 | 0,45              | violett/violet       | PM 21511              | 2   |
|   | 1.5          | 4-0 | 0,45              | ungefärbt/<br>undyed | PM 21512              | 2   |
|   | 2            | 3-0 | 0,45              | ungefärbt/<br>undyed | PM 21513              | 2   |

# CAPROLON RESORBA® P (LA/CL)

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

| Nadel / needle  | Stärke/gauge |     | Länge/length<br>m | Farbe/colour         | Best.-Nr.<br>code no. | Pckg.-Inhalt<br>box content<br>in Dtzd./dozen |
|---|--------------|-----|-------------------|----------------------|-----------------------|---|
|   | metric       | USP |                   |                      |                       |   |
| <br><b>DSM 18 F</b><br><br>PREMIUM-Schliff / PREMIUM cut        | 1.5          | 4-0 | 0,45              | violett/violet       | PM 2158               | 2   |
|   | 1.5          | 4-0 | 0,45              | ungefärbt/<br>undyed | PM 2159               | 2   |
| <br><b>DS 24</b>   | 1.5          | 4-0 | 0,45              | ungefärbt/<br>undyed | PM 2143               | 2   |
|   | 2            | 3-0 | 0,45              | ungefärbt/<br>undyed | PM 2144               | 2   |
|   | 3            | 2-0 | 0,45              | ungefärbt/<br>undyed | PM 2145               | 2   |
| <br><b>DSM 24</b><br>PREMIUM-Schliff / PREMIUM cut   | 1.5          | 4-0 | 0,45              | ungefärbt/<br>undyed | PM 21514              | 2   |
|   | 2            | 3-0 | 0,45              | ungefärbt/<br>undyed | PM 21515              | 2   |
|   | 3            | 2-0 | 0,45              | ungefärbt/<br>undyed | PM 21516              | 2   |
| <div style="background-color: #e0e0e0; padding: 10px;"> <p><b>3/8-Kreis schneidende Rundkörper (Diamond-Spitze), mit PREMIUM-Schliff = DRTM</b><br/> <i>3/8 circle, round bodied, trocar point, PREMIUM cut = DRTM</i></p> </div> |              |     |                   |                      |                       |   |
| <br><b>DRTM 10</b><br>PREMIUM-Schliff / PREMIUM cut  | 0.5          | 7-0 | 0,70              | violett/violet       | PM 2251               | 2   |
|   | 0.7          | 6-0 | 0,70              | violett/violet       | PM 2252               | 2   |
| <br><b>2 x DRTM 10</b><br>PREMIUM-Schliff / PREMIUM cut  | 0.5          | 7-0 | 0,70              | violett/violet       | PM 2253               | 2   |
|   | 0.7          | 6-0 | 0,70              | violett/violet       | PM 2254               | 2   |
| <br><b>DRTM 12</b><br>PREMIUM-Schliff / PREMIUM cut  | 0.7          | 6-0 | 0,70              | violett/violet       | PM 2255               | 2   |
|   | 1            | 5-0 | 0,70              | violett/violet       | PM 2256               | 2   |
| <br><b>2 x DRTM 12</b><br>PREMIUM-Schliff / PREMIUM cut  | 0.7          | 6-0 | 0,70              | violett/violet       | PM 2257               | 2   |
|   | 1            | 5-0 | 0,70              | violett/violet       | PM 2258               | 2   |