

NYLON monofil

Sterile Nadel-Faden-Kombination











Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
1/2-Kreis Rundkörper = HR F <i>1/2 circle, round bodied = HR F</i>						
● 2 x HR 6 F 	0.2	10-0	0,10	schwarz/black	5021	2
3/8-Kreis Rundkörper = DR, DR F <i>3/8 circle, round bodied = DR, DR F</i>						
● DR 4 	0.1	11-0	0,10	schwarz/black	50713	2
	0.2	10-0	0,10	schwarz/black	50714	2
● DR 5 F 	0.2	10-0	0,15	schwarz/black	5071	2
	0.3	9-0	0,15	schwarz/black	5072	2
● 2 x DR 5 F 	0.2	10-0	0,15	schwarz/black	5073	2
● DR 6 	0.2	10-0	0,15	schwarz/black	5074	2
	0.3	9-0	0,15	schwarz/black	5075	2
	0.4	8-0	0,15	schwarz/black	5076	2
● DR 10 F 	0.5	7-0	0,45	schwarz/black	5079	3
● DR 12 F 	0.7	6-0	0,75	schwarz/black	50711	3
	1	5-0	0,75	schwarz/black	50712	3
3/8-Kreis schneidender Rundkörper, mit PREMIUM-Schliff = DTRM <i>3/8 circle, trocar point, round bodied, PREMIUM cut = DTRM</i>						
⊙ DTRM 5  PREMIUM-Schliff / PREMIUM cut	0.2	10-0	0,15	schwarz/black	5251	2
1/2-Kreis schneidende = HS, mit PREMIUM-Schliff = HSM <i>1/2 circle, reverse cutting = HS, PREMIUM cut = HSM</i>						
▽ 2 x HSM 8  PREMIUM-Schliff / PREMIUM cut	0.3	9-0	0,45	schwarz/black	5121	2

NYLON monofil

Sterile Nadel-Faden-Kombination









Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
▼ HS 12 	1	5-0	0,45	schwarz/black	5111	3
▽ HSM 17  PREMIUM-Schliff / PREMIUM cut	1.5	4-0	0,45	schwarz/black	5126	3
3/8-Kreis schneidende = DS, DS F, mit PREMIUM-Schliff = DSM, DSM F <i>3/8 circle, reverse cutting = DS, DS F, PREMIUM cut = DS, DSM F</i>						
▽ DSM 6  PREMIUM-Schliff / PREMIUM cut	0.2	10-0	0,30	schwarz/black	5152	2
▽ DSM 7  PREMIUM-Schliff / PREMIUM cut	0.5	7-0	0,45	schwarz/black	5155	3
▽ DSM 11  PREMIUM-Schliff / PREMIUM cut	0.5 0.7	7-0 6-0	0,45 0,45	schwarz/black schwarz/black	5156 5157	3 3
▼ DS 13 	0.7 1 1.5	6-0 5-0 4-0	0,45 0,45 0,45	schwarz/black schwarz/black schwarz/black	5142 5143 5144	3 3 3
▽ DSM 13  PREMIUM-Schliff / PREMIUM cut	1 1.5	5-0 4-0	0,45 0,45	schwarz/black schwarz/black	51513 51514	3 3
▼ DS 16 	0.7 1	6-0 5-0	0,45 0,45	schwarz/black schwarz/black	5145 5146	3 3
▽ DSM 16  PREMIUM-Schliff / PREMIUM cut	1.5	4-0	0,45	schwarz/black	51517	3
▼ DS 18 	1 1.5 2 3 1.5 1 1.5 2	5-0 4-0 3-0 2-0 4-0 5-0 4-0 3-0	0,45 0,45 0,45 0,45 0,45 0,75 0,75 0,75	schwarz/black schwarz/black schwarz/black schwarz/black schwarz/black weiß/white schwarz/black schwarz/black schwarz/black	5147 5148 5149 51423 51414 51444 51445 51410	3 3 3 3 3 3 3 3

NYLON monofil

Sterile Nadel-Faden-Kombination












Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
▼ DS 18 F 	1.5	4-0	0,45	schwarz/black	51442	3
▽ DSM 18  PREMIUM-Schliff / PREMIUM cut	1.5 2 3 2	4-0 3-0 2-0 3-0	0,45 0,45 0,45 0,75	schwarz/black schwarz/black schwarz/black schwarz/black	51519 51520 51521 51522	3 3 3 3
▽ DSM 18 F  PREMIUM-Schliff / PREMIUM cut	1.5	4-0	0,45	schwarz/black	51539	3
▼ DS 21 	1.5 2	4-0 3-0	0,45 0,45	schwarz/black schwarz/black	51439 51440	3 3
▼ DS 24 	2 2	3-0 3-0	0,45 0,75	schwarz/black schwarz/black	51415 51416	3 3
▽ DSM 24  PREMIUM-Schliff / PREMIUM cut	2 3 2 3	3-0 2-0 3-0 2-0	0,45 0,45 0,75 0,75	schwarz/black schwarz/black schwarz/black schwarz/black	51523 51524 51535 51536	3 3 3 3
▼ DS 27 	3	2-0	0,45	schwarz/black	51419	3
▼ DS 30 	3.5 3	0 2-0	0,45 0,75	schwarz/black schwarz/black	51421 51454	3 3

NYLON monofil

Sterile Nadel-Faden-Kombination





Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
1/2-Kreis schneidende Spatula, PREMIUM-Schliff = HSPM, HSPM F <i>1/2 circle, spatulated, PREMIUM cut = HSPM, HSPM F</i>						
 HSPM 6 F  PREMIUM-Schliff / PREMIUM cut	0.2	10-0	0,15	schwarz/black	51914	2
	0.3	9-0	0,15	schwarz/black	51913	2
2 x HSPM 6 F  PREMIUM-Schliff / PREMIUM cut						
0.2	10-0	0,30	schwarz/black	5198	2	
2 x HSPM 7  PREMIUM-Schliff / PREMIUM cut						
0.2	10-0	0,30	schwarz/black	5193	2	
0.3	9-0	0,30	schwarz/black	5194	2	
3/8-Kreis schneidende Spatula, PREMIUM-Schliff = DSPM <i>3/8 circle, spatulated, PREMIUM cut= DSPM</i>						
 DSPM 4  PREMIUM-Schliff / PREMIUM cut	0.1	11-0	0,10	schwarz/black	52013	2
	0.2	10-0	0,10	schwarz/black	52014	2
DSPM 5  PREMIUM-Schliff / PREMIUM cut						
0.2	10-0	0,30	schwarz/black	5201	2	
2 x DSPM 5  PREMIUM-Schliff / PREMIUM cut						
0.2	10-0	0,30	schwarz/black	5202	2	
 DSPM 6  PREMIUM-Schliff / PREMIUM cut	0.2	10-0	0,15	schwarz/black	5206	2
	0.3	9-0	0,15	schwarz/black	52011	2
	0.2	10-0	0,30	schwarz/black	52015	2
2 x DSPM 6  PREMIUM-Schliff / PREMIUM cut						
0.2	10-0	0,30	schwarz/black	5208	2	

NYLON monofil

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
 2 x DSPM 7  PREMIUM-Schliff / PREMIUM cut	0.2	10-0	0,30	schwarz/black	5203	2
	0.3	9-0	0,30	schwarz/black	5204	2
<div style="background-color: #e0e0e0; padding: 5px;"> <p>1/4-Kreis schneidende Spatula, PREMIUM-Schliff = VSPM <i>1/4 circle, spatulated, PREMIUM cut = VSPM</i></p> </div>						
 2 x VSPM 7  PREMIUM-Schliff / PREMIUM cut	0.3	9-0	0,45	schwarz/black	5214	2

NYLON monofil

Einzelfadenpackungen, steril in vorgeschrittenen Längen

Sterile pre-cut suture, non-neededled

Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
metric	USP				
1.5	4-0	3 x 0,50	schwarz/black	X 335	2
2	3-0	3 x 0,50	schwarz/black	X 435	2
3	2-0	3 x 0,50	schwarz/black	X 535	2