











Resopren® PVDF-Faden

Sterile Nadel-Faden-Kombination











Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen	
	metric	USP					
1/2-Kreis Rundkörper = HR, HRS <i>1/2 circle, round bodied = HR, HRS</i>							
● HR 10 	0.7	6-0	0,45	blau/blue	77021	3	
● 2 x HR 10 	0.7	6-0	0,75	blau/blue	77022	3	
● 2 x HR 12 	0.7	6-0	0,75	blau/blue	77023	3	
	1	5-0	0,75	blau/blue	77024	3	
	1.5	4-0	0,90	blau/blue	77025	3	
● 2 x HR 15 	0.7	6-0	0,75	blau/blue	77026	3	
	1	5-0	0,75	blau/blue	77027	3	
	0.7	6-0	4 x 0,75	blau/blue	L 77028	1	
						Multi-L-Pack	
	1	5-0	4 x 0,75	blau/blue	L 77029	1	
					Multi-L-Pack		
● HR 17 	1	5-0	0,75	blau/blue	770210	3	
	1.5	4-0	0,75	blau/blue	770211	3	
	2	3-0	0,75	blau/blue	770212	3	
● 2 x HR 17 	1	5-0	0,90	blau/blue	770213	3	
	1.5	4-0	0,90	blau/blue	770214	3	
	2	3-0	0,90	blau/blue	770215	3	
● HR 22 	1.5	4-0	0,75	blau/blue	770216	3	
	2	3-0	0,75	blau/blue	770217	3	
● 2 x HR 22 	1.5	4-0	0,90	blau/blue	770218	3	
	2	3-0	0,90	blau/blue	770219	3	
● HR 27 	2	3-0	0,75	blau/blue	770220	3	
	3	2-0	0,75	blau/blue	770221	3	
● 2 x HR 27 	1.5	4-0	0,90	blau/blue	770222	3	
	2	3-0	0,90	blau/blue	770223	3	
	3	2-0	0,90	blau/blue	770224	3	

Resopren® PVDF-Faden

Sterile Nadel-Faden-Kombination











Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
● HRS 27 	3	2-0	0,75	blau/blue	770227	3
	3.5	0	0,75	blau/blue	770228	3
1/2-Kreis schneidende Rundkörper = HRT 1/2 circle, trocar point, round bodied = HRT						
▼ HRT 18 	1	5-0	0,75	blau/blue	77061	3
	1.5	4-0	0,75	blau/blue	77062	3
	2	3-0	0,75	blau/blue	77063	3
▼ 2 x HRT 18 	1.5	4-0	0,90	blau/blue	77065	3
	2	3-0	0,90	blau/blue	77066	3
	3	2-0	0,90	blau/blue	77067	3
▼ HRT 22 	1.5	4-0	0,75	blau/blue	77068	3
	2	3-0	0,75	blau/blue	77069	3
	3	2-0	0,75	blau/blue	770610	3
▼ 2 x HRT 22 	1.5	4-0	0,90	blau/blue	770611	3
▼ HRT 25 	2	3-0	0,75	blau/blue	770615	3
	3	2-0	0,75	blau/blue	770616	3
	3.5	0	0,75	blau/blue	770621	3
▼ 2 x HRT 25 	2	3-0	0,90	blau/blue	770617	3
	3	2-0	0,90	blau/blue	770618	3
3/8-Kreis Rundkörper = DR, DR F 3/8 circle, round bodied = DR, DR F						
● DR 10 	0.7	6-0	0,75	blau/blue	77074	3
● 2 x DR 10 	0.7	6-0	0,75	blau/blue	77075	3
	0.7	6-0	4 x 0,75	blau/blue	L 77076 Multi-L-Pack	1
● DR 10 F 	0.5	7-0	0,75	blau/blue	77071	3

Resopren® PVDF-Faden

Sterile Nadel-Faden-Kombination












Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
● 2 x DR 10 F 	0.5	7-0	0,75	blau/blue	77072	3
	0.5	7-0	4 x 0,75	blau/blue	L 77073 Multi-L-Pack	1
● DR 12 	1	5-0	0,75	blau/blue	770711	3
	1.5	4-0	0,75	blau/blue	770712	3
● 2 x DR 12 	1	5-0	0,75	blau/blue	770713	3
	1.5	4-0	0,75	blau/blue	770714	3
	1	5-0	0,90	blau/blue	770715	3
	1	5-0	4 x 0,90	blau/blue	L 770716 Multi-L-Pack	1
● DR 12 F 	0.7	6-0	0,75	blau/blue	77078	3
● 2 x DR 12 F 	0.7	6-0	0,75	blau/blue	77079	3
	0.7	6-0	4 x 0,75	blau/blue	L 770710 Multi-L-Pack	1
● DR 18 	1	5-0	0,75	blau/blue	770718	3
	1.5	4-0	0,75	blau/blue	770719	3
	2	3-0	0,75	blau/blue	770720	3
● 2 x DR 18 	1	5-0	0,90	blau/blue	770721	3
	1.5	4-0	0,90	blau/blue	770722	3
	2	3-0	0,90	blau/blue	770723	3
<p>3/8-Kreis schneidende Rundkörper = DRT, mit PREMIUM-Schliff = DRTM <i>3/8 circle, trocar point, round bodied = DRT, PREMIUM cut = DRTM</i></p>						
▽ DRTM 10  PREMIUM-Schliff / PREMIUM cut	0.5	7-0	0,75	blau/blue	77251	2
	0.7	6-0	0,75	blau/blue	77252	2
▽ 2 x DRTM 10  PREMIUM-Schliff / PREMIUM cut	0.5	7-0	0,75	blau/blue	77253	2
	0.7	6-0	0,75	blau/blue	77255	2
▽ DRT 12 	0.7	6-0	0,75	blau/blue	77081	3
	1	5-0	0,75	blau/blue	77082	3

Resopren® PVDF-Faden

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
 2 x DRT 12	0.7	6-0	0,75	blau/blue	77083	3
	1	5-0	0,75	blau/blue	77084	3
 DRTM 12	0.7	6-0	0,75	blau/blue	77257	2
	1	5-0	0,75	blau/blue	77258	2
 2 x DRTM 12 PREMIUM-Schliff / PREMIUM cut	0.7	6-0	0,75	blau/blue	772510	2
	1	5-0	0,75	blau/blue	772511	2
 DRT 18	1	5-0	0,75	blau/blue	77086	3
	1.5	4-0	0,75	blau/blue	77087	3
	2	3-0	0,75	blau/blue	77088	3
3/8-Kreis schneidende = DS, mit PREMIUM-Schliff = DSM, DSMF, DSMFX <i>3/8 circle, reverse cutting = DS, PREMIUM cut = DSM, DSMF, DSMFX</i>						
 DSM 7 PREMIUM-Schliff / PREMIUM cut	0.5	7-0	0,45	blau/blue	771518	3
	0.7	6-0	0,45	blau/blue	771519	3
 DSM 11 PREMIUM-Schliff / PREMIUM cut	0.7	6-0	0,45	blau/blue	771520	3
 DS 13	1	5-0	0,45	blau/blue	771411	3
 DSM 13 PREMIUM-Schliff / PREMIUM cut	0.7	6-0	0,45	blau/blue	771526	3
	1	5-0	0,45	blau/blue	771517	3
	1.5	4-0	0,45	blau/blue	771527	3
 DSMF 13 PREMIUM-Schliff "fein" / PREMIUM cut "fine"	0.7	6-0	0,45	blau/blue	77221	3
	1	5-0	0,45	blau/blue	77222	3
 DS 16 PREMIUM-Schliff / PREMIUM cut	1	5-0	0,45	blau/blue	771415	3
	1.5	4-0	0,45	blau/blue	771412	3
 DSM 16 PREMIUM-Schliff / PREMIUM cut	1	5-0	0,45	blau/blue	771525	3
	1.5	4-0	0,45	blau/blue	77151	3

Resopren® PVDF-Faden

Sterile Nadel-Faden-Kombination








Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen
	metric	USP				
 ▽ DSMF 16 PREMIUM-Schliff "fein" PREMIUM cut "fine"	1	5-0	0,45	blau/blue	77223	3
	1.5	4-0	0,45	blau/blue	77224	3
 ▼ DS 18	1	5-0	0,45	blau/blue	771401	3
	1.5	4-0	0,45	blau/blue	771402	3
	2	3-0	0,45	blau/blue	771403	3
	1.5	4-0	0,75	blau/blue	771416	3
	2	3-0	0,75	blau/blue	771417	3
	1.5	4-0	0,60	blau/blue	771404	2
	2	3-0	0,60	blau/blue	771405	2
				mit Intrakutanzubehör <i>Intracutaneous suture</i>		
				mit Intrakutanzubehör <i>Intracutaneous suture</i>		
 ▽ DSM 18 PREMIUM-Schliff / PREMIUM cut	1	5-0	0,45	blau/blue	77152	3
	1.5	4-0	0,45	blau/blue	77153	3
	2	3-0	0,45	blau/blue	77154	3
	1.5	4-0	0,75	blau/blue	771529	3
	2	3-0	0,75	blau/blue	77155	3
	3	2-0	0,75	blau/blue	771530	3
	1.5	4-0	0,60	blau/blue	77156	2
				mit Intrakutanzubehör <i>Intracutaneous suture</i>		
				mit Intrakutanzubehör <i>Intracutaneous suture</i>		
 ▽ DSM 18 F PREMIUM-Schliff / PREMIUM cut	1.5	4-0	0,75	blau/blue	771522	3
 ▽ DSMF 18 PREMIUM-Schliff "fein" PREMIUM cut "fine"	1	5-0	0,45	blau/blue	77225	3
	1.5	4-0	0,45	blau/blue	77226	3
	2	3-0	0,45	blau/blue	77227	3
	1.5	4-0	0,60	blau/blue	77228	2
	2	3-0	0,60	blau/blue	77229	2
				mit Intrakutanzubehör <i>Intracutaneous suture</i>		
				mit Intrakutanzubehör <i>Intracutaneous suture</i>		
 ▼ DS 21	2	3-0	0,45	blau/blue	771420	3
	3	2-0	0,45	blau/blue	771421	3
	2	3-0	0,75	blau/blue	771422	3
	3	2-0	0,75	blau/blue	771423	3
 ▽ DSM 21 PREMIUM-Schliff / PREMIUM cut	2	3-0	0,45	blau/blue	771534	3
	1.5	4-0	0,75	blau/blue	771523	3
	2	3-0	0,75	blau/blue	771524	3

Resopren® PVDF-Faden

Sterile Nadel-Faden-Kombination

Needle - suture - combination, sterile

Nadel / needle	Stärke Nr. gauge no.		Länge/length m	Farbe/colour	Best.-Nr. code no.	Pckg.-Inhalt box content in Dtzd./dozen	
	metric	USP					
▼ DS 24 	1.5	4-0	0,45	blau/blue	771418	3	
	2	3-0	0,45	blau/blue	771406	3	
	3	2-0	0,45	blau/blue	771419	3	
	2	3-0	0,75	blau/blue	771407	3	
	3	2-0	0,75	blau/blue	771408	3	
	2	3-0	0,60	blau/blue	771409	2	
					mit Intrakutanzubehör <i>Intracutaneous suture</i>		
	3	2-0	0,60	blau/blue	771410	2	
				mit Intrakutanzubehör <i>Intracutaneous suture</i>			
▽ DSM 24  PREMIUM-Schliff / PREMIUM cut	1.5	4-0	0,45	blau/blue	771531	3	
	2	3-0	0,45	blau/blue	77158	3	
	3	2-0	0,45	blau/blue	771532	3	
	2	3-0	0,75	blau/blue	77159	3	
	3	2-0	0,75	blau/blue	771510	3	
	2	3-0	0,60	blau/blue	771511	2	
					mit Intrakutanzubehör <i>Intracutaneous suture</i>		
	3	2-0	0,60	blau/blue	771512	2	
				mit Intrakutanzubehör <i>Intracutaneous suture</i>			
▽ DSMFX 24 	2	3-0	0,45	blau/blue	772210	3	
	2	3-0	0,60	blau/blue	772211	2	
					mit Intrakutanzubehör <i>Intracutaneous suture</i>		
▼ DS 27 	3	2-0	0,45	blau/blue	771425	3	
	3	2-0	0,75	blau/blue	771426	3	
▼ DS 30 	2	3-0	0,75	blau/blue	771424	3	
	3	2-0	0,75	blau/blue	771413	3	
	3.5	0	0,75	blau/blue	771414	3	
▽ DSM 30  PREMIUM-Schliff / PREMIUM cut	2	3-0	0,75	blau/blue	771513	3	
	3	2-0	0,75	blau/blue	771514	3	
	3.5	0	0,75	blau/blue	771515	3	
▼ DS 40 	3	2-0	0,75	blau/blue	771427	3	