




RESODURA[®] *matrix onlay*

Collagen Dura Substitute Membrane

RESODURA[®]*matrix onlay* is a conformable and resorbable membrane matrix engineered from highly purified type I collagen. It is indicated for use as an onlay (non-sutured) collagen membrane for the repair of the dura mater. RESODURA[®]*matrix onlay* requires only 30 seconds of hydration with sterile saline prior to implantation and can be cut in either the dry or hydrated state. It has a thickness similar to that of native dura, conforms to the contours of the brain and does not require suturing. Minimal tension sutures may be used if desired. It facilitates healing with balanced in vivo resorption of the implant and regeneration of host tissue. RESODURA[®]*matrix onlay* is supplied sterile in a double peel pouch and can be stored at room temperature.



- Highly purified type I collagen
- Conformable
- Implantable without sutures
- Leak resistant
- Controlled resorption for effective dural repair
- Excellent handling characteristics

Collagen Dura Substitute Membrane

DUM02525
2.5 x 2.5 cm

DUM02575
2.5 x 7.5 cm

DUM00505
5 x 5 cm

DUM07575
7.5 x 7.5 cm

DUM010125
10 x 12.5 cm

DUM0125175
12.5 x 17.5 cm

REF	SIZE	PACKAGE CONTENTS
DUM02525	2.5 x 2.5 cm	1 membrane
DUM02575	2.5 x 7.5 cm	1 membrane
DUM00505	5 x 5 cm	1 membrane
DUM07575	7.5 x 7.5 cm	1 membrane
DUM010125	10 x 12.5 cm	1 membrane
DUM0125175	12.5 x 17.5 cm	1 membrane



Collagen Matrix, Inc.

15 Thornton Road, Oakland, NJ 07436, USA

Tel: 201.405.1477, **Toll Free:** 888.405.1001

Fax: 201.405.1355

www.collagenmatrix.com